

How To Retube Your Ear Mould

How to change the tubing

** This may already be done for you.*



Pull old tubing out slowly and steadily.



Cut end off new tubing at a slant*...



...so you are left with a tapered end.



Push tapered end into earmould.



Thread tapered end through...



...until bend reaches earmould.



Cut close to earmould.



Line up with old tubing and cut...



...to same length.



Reattach hearing aid to earmould.



Correctly fitted tubing.

For video instruction on retubing visit:

www.pat.nhs.uk/hearing-aid-maintenance.htm

Your Click and Post Service

As a Pennine Care service user, we can now deliver hearing aid accessories straight to your door by post. This includes batteries and tubing.

Please visit:

<https://clickandpost.penninecare.nhs.uk/>

Troubleshooting Guide for Hearing aid Wearers

